

## iBioID Proteomics Protocol

Soderling Lab: Please reference *Science* **353**, 1123–1129 (2016)

Last updated: July 28<sup>th</sup>, 2016

**Notes:** Many of the below procedures should be done under “keratin free” conditions. This means use only powder-free Nitrile gloves (NOT Latex) in a laminar flow hood.

### Specialized Reagents to order (name: company, catalog #)

Biotin: Sigma; B4501

Ammonium Bicarbonate (Ambic): Sigma; A6141

Deoxycholate: Fisher; BP349-100

Protease Inhibitors: RPI; (AEBSF-A20010)(Leupeptin-L22035)(PepstatinA-P30100)

HA antibody resin: Sigma; A2095

Low protein binding 1.5ml eppendorf tubes: Sigma; Z666505

NeutrAvidin: Thermo; 29201

Slide-A-Lyzer MINI Dialysis 3,500 MWCO: Thermo; 69550

RapiGest SF Surfactant: Waters; 186001861

Spin Cups-Cellulose Acetate Filter: Thermo; 69702

pyruvate carboxylase Ab: Invitrogen; A21970

MAP2 Ab: GeneTex; GTX49163

Protein G agarose (capacity 20mg IgG/mL): Millipore; 16-266

Mass spec compatible water: Honeywell; LC365-1

BCA protein quantification kit: Thermo; 23225

### Equipment needed:

Micro Tissue Grinder or dounce homogenizer (2ml)

Sonicator

Ultracentrifuge

### Solutions:

- Make 1M Ammonium Bicarbonate (Ambic) in keratin free condition.
- Make 2%SDS, 10mL in keratin free condition.
- Make 1M NaCl, 10mL in keratin free condition.
- Make 1% TritonX100, 1% deoxycholate, 25mM LiCl, 10mL in keratin free condition.
- Make 50mM Ambic in keratin free condition.
- Make 0.1% RapiGest / 1mM biotin / 50mM Ambic (200uL of 5mM biotin/50mM Ambic solution, 40uL of 1M Ambic, 760uL of dw in 1 vial RapiGest) in keratin free condition.
- Make 50mM Ambic, 500mL x 2. Keep the beaker in the cold room a day before the experiment.
  - [Ambic FW=79.06, to make 50mM 1.976g/500mL]

- Lysis(-RIPA) buffer [50mM Tris/Cl pH7.5, 150mM NaCl, 1mM EDTA]
- 2x RIPA buffer [50mM Tris/Cl pH7.5, 150mM NaCl, 1mM EDTA, 0.4%SDS, 2%TritonX100, 2% deoxycolate]

### **Day 1**

#### Preparations:

- 1) Per sample: Pre-wash 20uL NeutrAvidin slurry with 2xRIPA buffer, 5 times
- 2) Per sample: Incubate 20uL of pyruvate carboxylase Ab with 10uL of Protein G agarose in 150uL of 0.1% TritonX-100/PBS + PIs, o/n at 4 degrees with rocking.
- 3) Per sample: Incubate 20uL of MAP2 Ab with 10uL of Protein G agarose in 150uL of 0.1% TritonX-100/PBS + PIs, overnight at 4\* with rocking.

#### Procedures:

- 1) Set aside 3 eppendorf tubes (1.7 mL) per sample
- 2) Take tissue samples from biotin labeled mice from liquid nitrogen and keep on dry ice
- 3) Homogenize sample with 300 uL of Lysis-R buffer (including protease inhibitors). This should be done using 15 strokes with a glass dounce homogenizer (2ml, Pierce).
- 4) Transfer to an eppendorf tube, add 350 uL of 2x RIPA buffer and spin in microfuge at 15,000xg for 10 min at 4°C.
- 5) Sonicate sample at level 8, 7 sec.; Cool on ice for 10 min, Repeat Sonication and ice for a total of 3x.
- 6) Centrifuge at 15,000 rpm for 10 min, Transfer the supernatant to a new eppendorf tube.
- 7) Centrifuge in the ultracentrifuge (TLA-55, 40,000 rpm [100,000xg], 30 min, 4°C).
- 8) Transfer supernatant to a new eppendorf tube and add SDS to a final 1% concentration.
- 9) Briefly vortex and spin @ RT in microcentrifuge at 15,000xg for 10min.
- 10) Boil (95°C for 5 min), Cool on ice
- 11) Microcentrifuge at 15,000xg for 30 min at 4°C
- 12) Transfer supernatant to a 15 mL conical tube and add protease inhibitors
- 13) Incubate with 17 uL of washed NeutrAvidin beads overnight at 4\* with rocking.

### **Day 2**

#### Preparations:

- 1) Pre wash the HA beads (100uL beads) in keratin free condition in low binding tube. Wash with 1x RIPA buffer first and wash with 50mM Ambic x 5 in keratin free condition.
- 2) Pre wash the spin cups (your samples x2) in 50mM Ambic x 5 keratin free condition.
- 3) Pre wash 1 dialysis unit per sample (Slide-A-Lyzer MINI Dialysis 3500)

MWCO). Dialyze against 500mL of 50mM Ambic a few hours. Check no leak. Rinse 3 times with 50mM Ambic.

**4)** Wash carboxylase Ab and MAP2 Ab on Protein G beads from yesterday with 1x RIPA buffer first and wash with 50mM Ambic x 5 in keratin free condition.

**5)** Pre-heat PCR tubes with RapiGest solution (60°C, 2hrs) and rinse with 50mM Ambic x 5 keratin free condition.

#### Procedures:

**1)** Discard supernatant from overnight NeutrAvidin pulldown, and Transfer beads to a low-protein binding tube.

**2)** Wash beads with 2% SDS x 2 (10 min each) at room temperature.

**3)** Wash beads with 1% TritonX100, 1% deoxycholate, 25 mM LiCl x2 (10 min each) at RT, Wash with 1 M NaCl x 2 (10 min each) at RT (500 uL each wash).

**4)** Wash with 50 mM Ambic x 5 (10 min each) @ RT (500 uL).

**5)** Rinse once with cold 0.1% RapiGest/2mM biotin/ 50 mM Ambic and transfer to PCR tubes prepared above.

**6)** Spin down in microcentrifuge, discard supernatant, and add **40 uL** of 0.1% RapiGest/ 2mM biotin/50 mM Ambic.

**7)** Keep at 60°C for 2 hours.

**8)** Transfer to spin column (Pierce Spin Cups-Cellulose Acetate Filter) and spin down for 30 s at 3381xg in a low binding tube.

**9)** Transfer the flow-through to a low binding tube with 20 uL of HA beads and 6 uL of PC/MAP2 Ab conjugated beads.

**10)** Incubate overnight at 4°C (while shaking on nutator).

### Day 3

#### Procedures:

**1)** Transfer sample to a spin column (Pierce Spin Cups-Cellulose Acetate Filter) with low binding tube and spin down for 30 sec at 3381xg.

**2)** Transfer the supernatant to a "Slide A lyzer MINI" and dialyze against 500 mL of 50mM Ambic twice (3 hours each).

**3)** Transfer to a low-protein binding tube.

**4)** Measure the protein concentration using BCA kit from Pierce.

**5)** Eluted protein is next lyophilized to dryness and resuspended in 20 uL of 1X SDS loading buffer supplemented with 10 mM dithiothreitol.

**6)** All 20 uL is loaded onto an Invitrogen NuPAGE 4-12% SDS-PAGE gel and run for approximately 5 min to electrophorese all proteins into the gel matrix. The entire MW range was then excised in a single gel-band and subjected to standardized in-gel reduction, alkylation, and tryptic digestion.

**7)** Proceed to Mass Spectrometry.