


# Use of Serial Multi-Template Liquid Biopsies in Triple Negative Breast Cancer Monitoring

Paul Y. Song, Jill Simmons, Paul W. Dempsey, Erich Klem, Chris Carter, Shahram Tahvilian, Brian Scollick, William Strauss  
Cynvenio Biosystems Inc., Westlake Village, CA

## Background

- Triple negative breast cancer (TNBC) remains the most aggressive sub-type of breast cancer in which up to 1/3 of all patients will relapse early and distantly.
- There remains no proven method to monitor and detect early recurrence.
- This ongoing study examines the use of serial longitudinal liquid biopsies using a multi-template approach that identifies and analyzes cell-free DNA (cfDNA) and circulating tumor cells (CTCs) as a solution.



**Figure 1. Multi-template approach via ClearID liquid biopsy**

ClearID was designed to analyze both circulating tumor cells (CTCs) and cfDNA alongside the patient's healthy DNA (germline) to provide a clear distinction between tumor and germline.

## Results

163 patients ( <b>77.25%</b> )	had evidence of mutations in at least one sample
48 patients ( <b>24.3%</b> )	had <b>no</b> evidence of mutations in any sample
59 patients ( <b>30.10%</b> )	had mutations on <b>two consecutive</b> blood draws
18 patients ( <b>9.18%</b> )	had the <b>same</b> mutation on <b>two consecutive</b> blood draws
1 patient ( <b>0.51%</b> )	had a mutation on <b>both</b> CTC and cfDNA on two consecutive draws

### Results: Recurrences

Recurrences thus far	<b>11 patients (5.68%)</b>
Types of recurrences	10 breast, 1 ovarian
Average time from treatment to time of recurrence	<b>21 months</b>


**All recurrences to date have had at least two consecutive draws with mutations in either cfDNA, CTCs, or both.**

<b>36%</b>	have had the same mutation on consecutive draws
<b>18%</b>	have had the same mutation on consecutive draws in both cfDNA and CTCs.


<b>Total Number of Samples</b>	<b>826</b>
% of cfDNA samples with mutations	22.84%
% of CTC samples with mutations	19.84%
% of samples with mutations on both cfDNA and CTCs	4.28%

**Table 3. Sample collection characteristics**

**TNBC Mutation Frequencies as of May 15th, 2017 (excluding mutations expressing p.C1003\* and p.L56S)**



**TP53 Mutation Expressed Protein Frequencies as of May 15th, 2017**



**Fig 2 and 2a. Most commonly seen mutations**

## Methods

- **211** patients with confirmed diagnosis of triple negative breast cancer and within 3 years of completion of therapy and in full remission were enrolled. **ClearID** liquid biopsies were performed every three months to look for cancer specific genomic mutations not seen in germline controls, but in cfDNA and CTCs using a custom 27-gene breast cancer panel.
- Reporting threshold was 1.0%.
- Prior validation of our gene panel demonstrated a false positive rate of 0.001 - 0.0007% in normal cfDNA and cell-based controls.

AKT1	3	BRAF	3	EGFR	12	GATA3	16	KRAS	4	MAP2K4	15	NCOR1	102	PTEN	15	RYR2	47
ARID1A	6	CDH1	45	ERBB2	7	HRAS	2	LRP1	68	MAP3K1	19	NRAS	2	RB1	35	RYR3	20
ATM	4	CDKN2A	6	ESR1	27	KMT2C	11	LRP2	101	MED12	7	PIK3CA	14	RUNX1	11	TP53	19

**Table 1. Custom 27-gene panel showing number of amplicons per gene**

- 621 amplicons targeting 27 genes
- FFPE, CTC and cfDNA compatible
- Average size: 134bp
- Coverage of 52.5kb
- Range size: 107-202bp

<b>Patients (N)</b>	<b>211</b>	
Median Age (Years)	53 (28 -84)	
Average time from Treatment	25 mos.	
Stage	I	35.7%
	II	45.2%
	III	16.6%
	IV	0.01%
	NA	0.02%
Average # of blood draws	4	

**Table 2. Patient demographics**

## Conclusions

- **The majority (77.25%) of TNBC patients in our study post-treatment were found to have genomic mutations in at least one liquid biopsy sample - either cfDNA, CTC, or both.**
- **Of those patients found to have a mutation on any one sample, only 6.7% developed a recurrence and 93.3% remain NED.**
- **Unlike patients who remained NED, patients who recurred displayed persistent evidence of mutations in cfDNA, CTCs, or both in serial samples.**
- **When patients have any mutations on two successive draws, they are 2.7x's more likely to recur compared to those with only one positive sample.**
- **When patients have the same mutation on two successive samples, they are 3.3x's more likely to recur compared to those with only one positive sample.**
- **This study highlights the need for serial longitudinal monitoring using a multi-template approach to better fully understand patient specific tumor biology and kinetics.**
- **This study remains open as we anticipate several more recurrences in the coming months in patients currently exhibiting the same mutation in two consecutive samples who remain NED at this time.**