

Labeling Kinesin with Benzylguanaine (BG)-oligonucleotide

Here, microtubules are used to affinity purify kinesin labeled with benzylguanaine (BG)-oligonucleotides from excess label. This also results in the removal of dead kinesin motors.

Preparing Microtubules (MTs):

- 1) Thaw 20ul aliquot of tubulin (20 mg/ml) at 37°C.
- 2) Immediately place on ice.
- 3) Polymerize tubulin by mixing together:
 - 20 ul 13 mg/ml tubulin
 - 20 ul 2X polymerization mix
- 4) Incubate at 37°C for 20 min.

From this point on, the MTs should never be placed ice as microtubules depolymerize in the cold.

- 5) Taxol stabilize the MTs. The final concentration of taxol should be equimolar to the tubulin concentration.
 - Taxol binds to tubulin in a 1:1 complex. Taxol (sold by Sigma as paclitaxel) is made to a 10mM stock in DMSO.
 - Adding taxol directly before the MTs have polymerized causes tubulin to precipitate (see Mitchison lab website for more details: <http://mitchison.med.harvard.edu/protocols.html>)
- 6) Mix together:
 - 40 ul MTs polymerization mix from previous step
 - 40 ul 1X BRB80 containing the appropriate amount of taxol.
- 7) Incubate at 37°C for at least 10 minutes.

Store this at room temperature! All of the following steps should be done at room temperature and contain taxol.

- 8) Separate tubulin dimer from polymerized MTs. Pipette MTs from step 7 onto a 60% glycerol (or 40% sucrose) cushion
 - The glycerol cushion should be an = or great volume than the MT volume.
- 9) Spin 50,000 rpm for 15 minutes at room temp in a TLA100.
- 10) Carefully remove glycerol cushion and then resuspend MT pellet in BRB80 (or desired dynein buffer) + 20uM taxol, 1mM DTT.
- 11) Measure OD280 to determine MT concentration (5ul MTs into 95ul 6M guanidine HCl).

Buffers

2X Polymerization Mix:

- 100 ul 5X BRB80
- 0.5 ul 1M DTT
- 50 ul 100% DMSO
- 5 ul 100 mM GTP
- 5 ul 100 mM MgCl₂
- 89.5 ul filter sterilized ddH₂O

Glycerol Cushion

- 60% glycerol
- 1X BRB12
- 1mM DTT
- 20uM Taxol

BRB80 Buffer

- 1X BRB80
- 1mM DTT
- 20uM Taxol

BRB12 Buffer

- 1X BRB12
- 1mM DTT
- 20uM Taxol

5X BRB80

- 400mM PIPES (pH6.8)
- 10mM MgCl₂
- 5mM EGTA

5X BRB12

- 60mM PIPES (pH 6.8)
- 10mM MgCl₂
- 5mM EGTA

Labeling the motor by binding and releasing from microtubules:

- 1) Mix the following adjusting for molar ratios accordingly (molar ratios):
 - 50ul of 3.5uM K560 (1x)
 - 22ul of ~240uM BG-oligo (20-30x)
- 2) Incubate at room temperature then add the following in this order:
 - 0.85ul of 0.1M AMPPNP (1mM)
 - 0.17ul of 10mM taxol (20uM)
 - 0.85ul of 100x Apyrase (1x ~ 6.6 units/ml)
 - 11.5 ul of 50uM MTs (3-4x) (Taxol should be in the buffer before MTs are added).
- 3) Volume is ~85ul making final concentrations of all the components
 - ~2uM K560.
 - ~60uM BG-oligo
 - ~1mM AMPPNP
 - ~20uM Taxol
 - ~6.6 units/ml Apyrase
 - ~6.7uM MTs
 - i) The more concentrated the components are, the more efficient the labeling reaction.
- 4) Incubate at room temperature for 15 min after mixing the solution by flicking the tube gently.
- 5) Centrifuge through a cushion in TLA100, 50krpm, 25C, 15min in a TLA 100.
- 6) Remove the top layer (S1) carefully taking a small amount of cushion.
 - Take 4ul of S1 for a gel sample.
- 7) Wash the tube and above the cushion by adding BRB12 to the top and removing.

Repeat 3X total removing a bit of cushion each time but leaving the pellet completely covered in cushion.

- 8) Remove the cushion carefully without disturbing pellet.
- 9) Wash the pellet 3X with BRB12 avoiding disturbing the pellet.
- 10) Resuspend the pellet in 50ul of Release buffer
 - take 3ul for a gel as P1.
- 11) Incubate release reaction for 5 min at room temperature.
- 12) Pellet MTs at 80krpm, 6 min, 25C with the TLA 100.
- 13) Remove supernatant and add sucrose to a final concentration of 10%
 - take 3ul for a gel as S2.
- 14) The supernatant supplemented with sucrose can be aliquoted, flash frozen and stored at -80C. This contains the purified kinesin motor labeled with oligo with minimal free BG-oligo. There will be a mixture of labeled and unlabeled kinesin if the reaction has not gone to completion.
- 15) Resuspend the pellet in 50ul of buffer and take a sample for P2.
- 16) Run gel of all samples with an appropriate standard to measure labeling efficiency (gel shift from ~8kda oligo) and yield (i.e. final concentration of S2 sample).

Buffers

Release Buffer

- 30ul 2M KCl
- 40ul 5x BRB80
- 0.2ul 1M DTT
- 10ul 100mM ATP*MgCl₂
- 0.4ul 10mM Taxol
- 119.4 ul ddH₂O

BRB12 Cushion

- 60% Glycerol
- 1x BRB12
- 1mM DTT
- 20uM Taxol

BRB12 Buffer

- 1xBRB12
- 1mM DTT
- 20uM Taxol