

StimRouter[®]

Peripheral Nerve Stimulation System

For the Treatment of Chronic Pain



Agenda



Bioness Introduction




Technical Overview



Patient Outcomes



Clinical Study



Referral Pathway



Next Steps

Mission Statement

“**Bioness** provides neuromodulation technologies that **help improve lives** and restore function for those living with neurological deficit and peripheral pain.”

Born from the Alfred Mann Foundation



- ✦ Cochlear Implants - Advanced Bionics
- ✦ Retinal Implants - Second Sight
- ✦ Insulin Delivery - Medtronic
- ✦ Cardiac Rhythm Management - St. Jude
- ✦ Functional Rehabilitation - Bioness

Technologies designed to allow:

“The blind to see. The deaf hear. The lame walk.”

Alfred Mann

Bioness Portfolio -A Brand You Know

L300 FES System
Foot Drop System
2006



L300 Plus FES System for
Thigh Weakness
2010



H200 Hand
Rehabilitation System
2004



Vector
Gait and Safety System
2013



L300 Go FES System
2017



Bioness Integrated Therapy
System (BITS)
2014



StimRouter Neuromodulation
Therapy for Pain
2015



❄ 7 FDA-cleared medical devices

❄ Award Winning Technologies

❄ 170 Licensed Patents

StimRouter®

Chronic Peripheral Nerve Pain Management



Common Clinical Applications



Failed surgery pain –
knee, hip, back



CRPS



Nerve compression,
injury or trauma



Back
pain



Post-stroke shoulder
pain



Spinal cord injury
pain



Foot/neuroma
pain



Post-amputation
pain

Treatment Options | Chronic Peripheral Nerve Pain



Rx Meds/
Opioid Use



Physical
Therapy



Nerve Block/
Ablation



Spinal Cord
Stimulation



Do
Nothing

Neuromodulation



“The alteration of nerve activity through targeted delivery of a stimulus, such as electrical stimulation or chemical agents, to specific neurological sites in the body.”

ins
International
Neuromodulation
Society

 Bioness

StimRouter®

Pain Management: Spinal Cord Stimulators

Completely implanted system

Lead placed in dorsal column of spine tunneled to an implantable pulse generator

Rechargeable & primary cell requires surgery to replace battery



Spinal Cord Stimulators vs PNS



70%

Number of Neuromodulation Procedures treated with SCS

6K

Off-Label SCS procedures used to treat peripheral nerve pain per year

**Diverse etiologies and pathologies;
Many implant locations not easily
accessible by SCS**

Post Stroke Shoulder Pain

>10M

Number of Stroke Survivors
in the United States

800K

Number of Strokes Each
Year in the United States

Up to
~60%

Number of Stroke Patients
Reporting Pain each year

Stroke Recovery Challenges:

PSSP is difficult to treat and limits recovery potential.
Early intervention is key.



Post Stroke Shoulder Pain

Shoulder pain can occur as early as 2 weeks post-stroke but typically occurs within 2-3 months post-stroke.

Contributing Factors

- ✿ Subluxation
- ✿ Contractures
- ✿ Rotator cuff injury
- ✿ Spastic muscle imbalance of the glenohumeral joint

Limits motor and functional recovery, ADLs. Decreases QoL.



StimRouter[®]

**Target Pain
at its Source**






CHRONIC PAIN



StimRouter

 **Bioness**

StimRouter Neuromodulation System



FDA cleared, minimally-invasive implant designed to treat **chronic pain of peripheral nerve origin**, below the cranial facial region. The minimally-invasive lead implant procedure is performed **under local anesthesia** through a small incision.

Powered externally through the skin to stimulate the peripheral nerve with a small, focal electrical field - interrupting the pain signal to alleviate pain.

Puts **patients in control** of their pain with a handheld, wireless Patient Programmer.

System Components



External Pulse Transmitter (EPT) is removed when not in use & recharged nightly;
Gel electrode patch replaced every 2-4 days

Patient Programmer controls stimulation and commands EPT to run up to 8 customized stimulation programs

Lead contains 3 stimulation electrodes, integrated receiver and anchor; 15cm long, fully implanted

StimRouter[®]

© 2018

StimRouter Lead



- ✿ Only component of the system that is implanted
- ✿ Flexible, durable 15cm lead with integrated receiver
- ✿ Anchoring mechanism designed to prevent migration
- ✿ Three, single channel electrode contacts



External Pulse Transmitter (EPT)



- ✱ Transmits electrical field stimulation to receiving end of lead that is implanted under the skin
- ✱ Programmable: Stores up to 8 stimulation programs
- ✱ Rechargeable: Can operate ~2 days on single charge
- ✱ Attaches to disposable Electrode Patch
- ✱ Gel patch adheres to skin to properly position EPT

StimRouter®

© 2018

Patient Programmer



- ✿ Wirelessly controls the EPT
- ✿ Turns on/off; Adjusts intensity +/-
- ✿ Allows patients to control stimulation, manage programs & intensity
- ✿ Tracks compliance & usage
- ✿ Visual and auditory indicators
- ✿ Rechargeable battery

StimRouter®

© 2018

Contradictions/Precautions

Precautions

Patients who have a cancerous lesion near the target stimulation point or near the electrode patch. Patients exposed to diathermy. Patient exposed to therapeutic ultrasound.

Exclusions

Demand Cardiac Pacemaker or Defibrillator. Metallic implant in the immediate area of the implant (e.g. shoulder). Patients with Central Pain Syndromes.

A 3D medical animation showing a robotic hand holding a dark grey surface probe. The probe is being used to touch the skin of a human arm, specifically the area of the axilla. The background is a teal color with a white surgical drape on the left. The text 'Implant Procedure Animation' is overlaid in white on the top left.

Implant Procedure Animation

A surface probe can be used to locate motor point of the Axillary Nerve.

Clinical Study Results

Primary Efficacy

At 3-months the group receiving StimRouter treatment demonstrated a **statistically significant improvement** in pain as compared to the control group ($p < 0.0001$).

Primary Safety

No serious adverse events related to the device were reported during the duration of the study (12-months).

Secondary Outcomes

The treatment group had significantly more **favorable outcomes related to quality of life** and satisfaction as compared to those in the control group.

StimRouter®

© 2018

Clinical Study Results

50%

of the treatment group rated their **satisfaction 8 or higher** on a 10 point scale

53%

of the treatment group rated their global impression of change in activity limitations, symptoms, emotions and **overall quality of life** related to their painful condition between 5-7 on a 7-point scale

31%

of the treatment group rated their **satisfaction at a 10** on a 10-point scale

Targeted Peripheral Nerves



Arm

Ulnar

Median

Radial

Axillary

Suprascapular

Trunk

Ilioinguinal

Intercostal

Pudendal

Iliohypogastric

Coccygea

Genitofemoral

Superior Cluneal

Leg

Saphenous

Tibial



Femoral

Femoral Cutaneous

Sural

Bioness Patient Survey

15 different nerve locations | Carry-over effect, decreases utilization over time | Data up to 5 years post-clinical study



Old Pain Algorithm

1

Temporary
pain relief

2

Temporary,
inconsistent
outcomes

3

Sometimes painful,
not durable, more
damage to nerve

4

Invasive surgery,
not reversible

5

Invasive surgery,
difficult for many
peripheral nerve
locations



New Pain Algorithm

1

Diagnostic Block

2

PNS/StimRouter

3

Periph. Nerve
Surgery

4

Cryo/Thermal
Ablation

5

SCS/ITP



Why StimRouter?

Target Peripheral Nerve Pain at its Source

20+ locations implanted, 78% responder rate, 68% pain relief overall, reduced use of Rx medications

Minimally-Invasive

Two small incisions, often a 30 minute procedure, local anesthesia

Durable Pain Relief

Compared to injections and other temporary interventions, StimRouter may provide lasting relief

Reversible and MRI Conditional

MRI scans permitted under certain conditions; Single incision explant procedure if required

Patient Testimonials



LUCINA: “When the doctor turned it on, **the pain went away** and now all I feel is a slight tingling. Now I can **move my arm, lift it, and do exercises** with the help of my therapist and the pain is totally gone.”

West Palm Beach, FL




KAYE: “I was surprised as to how well the StimRouter worked. **It was a life-saver.** When you’re an active person like I am and all of a sudden your body doesn’t work and you’re in constant pain, it’s very tough. Now, I don’t have pain anymore.”

Redondo Beach, CA

Clinical Centers of Excellence

Stanford University
Mount Sinai
UCLA
University of UT
VCU
Cleveland Clinic
Johns Hopkins



StimRouter®

Deliver StimRouter Outcomes to Your Patients

Bioness Support for Implanting Physicians



Training format
customized to
your needs



Assistance with
technology
onboarding/VAC
process



Onsite support for
first implants and
patient programming
appointments



Ongoing
reimbursement
support
(documentation,
pre-auths, appeals)



Marketing support
for patient referrals

Deliver StimRouter Outcomes to Your Patients

Role of the Referring Physician



- 1 Refer to Bioness for pre-screening and connection to a local implanter; Form 1 authorizes Bioness to contact patient
-Or refer patients to www.stimrouter.com/dtcinquiries
- 2 Implanting physician will evaluate for PNS treatment
- 3 Patient referred back to you post-implant
- 4 Bioness to track user outcomes



Kaye Gagnon



Bioness



StimRouter®

Referenced Material

Evidence-Based Review of Stroke Rehabilitation. Chapter 11: Painful Hemiplegic Shoulder.

http://www.ebrsr.com/sites/default/files/chapter11_hemiplegicshoulder_final_16ed.pdf

Deer T, et al. 2016. *Neuromodulation*. 19:91-100.

Poduri KR. Shoulder pain in stroke patients and its effect on rehabilitation. *J Stroke Cerebrovascular Dis*. 1993. 3:261-6.

Coskun Benlidayi I, Basaran S. Hemiplegic shoulder pain: a common clinical consequence of stroke. *Pract Neurol*. 2013 Aug 12. [\[Medline\]](#).

Adey-Wakeling Z, Arima H, Crotty M, et al. Incidence and associations of hemiplegic shoulder pain poststroke: prospective population-based study. *Arch Phys Med Rehabil*. 2015 Feb. 96 (2):241-247.e1. [\[Medline\]](#).

Teasell RW. The painful hemiplegic shoulder. *Physical Medicine and Rehabilitation: State of the Art Reviews*. 1998. 12(3):489-500.

Forster A. painful hemiplegic shoulder: physiotherapy treatment. *Rev Clin Gerontol*. 1994. 4:343-8.

Najenson T, Yacubovich E, Pikielni SS. Rotator cuff injury in shoulder joints of hemiplegic patients. *Scand J Rehabil Med*. 1971. 3(3):131-7. [\[Medline\]](#).

Shoulder Pain After Stroke: A Prospective Population-Based Study , 2006 <http://stroke.ahajournals.org/content/38/2/343.full>

Spinal Cord Stimulation's Role in Managing Chronic Disease Symptoms, INS 2013 <http://www.neuromodulation.com/spinal-cord-stimulation>

Van Ouwenaller, C, Laplace, P. M., & Chantraine A. (1986) Painful shoulder in hemiplegia. *Archives of Physical Medicine and Rehabilitation*, 67, 23–26.

StimRouter®

Thank you.

Questions?

1.800.211.9136

info@stimrouter.com

www.stimrouter.com

