

BrainXell Astrocyte Monoculture

Protocol (v10.0)

CONTENTS

- One vial cryopreserved human astrocytes (BX-0600-XX-XX)
- BrainFast Astro (BrainXell #BX-2600) formerly known as Astrocyte Supplement
- BrainFast SK (BrainXell # BX-2020) formerly known as Supplement K

STORAGE

- Immediately transfer the cryovial of astrocytes to a liquid or vapor nitrogen storage system.
- BrainFast supplements should be stored at -20°C (6 months) or -80°C (18 months). Return vial to -20°C between each time of use to maintain stability.

ADDITIONAL MATERIALS NEEDED

- DMEM/F-12 Medium, +L-glutamine +HEPES (Thermo Fisher Scientific #11330032)
- Neurobasal Medium (Thermo Fisher Scientific #21103049)
- N-2 Supplement (Thermo Fisher Scientific #17502048)
- GlutaMAX Supplement (Thermo Fisher Scientific #35050061)
- Fetal Bovine Serum (Corning #35016CV)
- PDL-Coated 96-Well Plates (Refer to <u>BrainXell Culture Plates PDL Coating Protocol v10.0</u>)

PROCEDURE

Day O Seeding Preparation

To seed (1) full 96-well plate you will need 2.2 - 2.8 million live cells. Additional cells may be used depending on individual research needs. Review the Certificate of Analysis (CoA) for viability and seeding information.

- Gather the components for the Basal Medium according to the recipe in the table below (see Media Compositions section below).
- 2. Working in a cell culture hood (biological safety cabinet), combine all Basal Medium components in an appropriately sized sterile container. Allow the Basal Medium to equilibrate to room temperature for at least 15 minutes. Do not warm the medium in a 37°C water bath. Culture plates should also be at room temperature prior to use.
- 3. Prepare one 50-mL conical tube: add 3 mL of Basal Medium.
- 4. Prepare one microcentrifuge tube: add 25 µL of Trypan Blue solution for cell counting.

Day 0 Seeding the Astrocytes

- 5. Remove the cryovial from the liquid nitrogen and place in a 37°C water bath. To minimize contamination, avoid submerging the cap. Gently move the vial within the bath to increase the rate of thawing.
- 6. As soon as the last of the ice melts, which will take \sim 75-90 seconds, remove the cryovial from the water bath. Disinfect the vial by spraying it with 70% ethanol before transferring it into the cell culture hood.

- 7. Using a P1000 pipette, transfer 500 μ L Basal Medium from the prepared 50-mL conical tube to the cell cryovial at a rate of ~2-3 drops/sec. This process should take about 10 seconds.
- 8. Gently transfer all contents from the cryovial (~1 mL total) back to the same 50-mL conical tube.
- 9. Centrifuge cells at 1700 rpm (465xg) for 5 mins.
- 10. Carefully aspirate supernatant, leaving the cell pellet, and resuspend cells in 950 μ L fresh Basal Medium. Mix thoroughly by gently pipetting up and down using a P1000 pipette.
- 11. Based on the CoA value (Viable Cells/Vial), resuspend the cells to a concentration of 1.0×10^6 live cells/mL by slowly adding additional Basal Medium to the existing ~ 1 mL in the tube.
 - a. Example: 2.3 million Viable Cells/Vial is diluted to 2.3 mL total volume.
- 12. Count the cells: gently swirl the conical tube and pipette up and down 3-5 times to ensure cells are evenly suspended in medium. Transfer 25 μ L of cell suspension to the microcentrifuge tube prefilled with 25 μ L Trypan Blue solution from Step 4 and pipette up and down a few times to mix. Count the number of viable and dead cells using a hemocytometer. Determine the live cell concentration (live cells/mL) and viability.
- 13. Calculate volumes needed to make the Astrocyte Seeding Suspension. A typical seeding density is 20,000 -25,000 viable cells/100 μ L/well for a 96-well plate (\sim 62,500 78,000 viable cells/cm²). Recommended seeding density may vary based on lot number; refer to the CoA for lot-specific seeding information. Dilute the cells to the desired seeding concentration based on the Trypan Blue cell count.

Example of dilution calculations.

Actual Live Cell Concentration	Target Seeding Density	Total Seeding Volume Needed	Volume of Cell Suspension Volume	Volume of Basal Medium
Ex: 1.1x10 ⁶ viable cells/mL (Step 12)	20,000 viable cells per 100 μ L/well = 200,000 cells/mL	11 mL (Seeding 1 plate)	Mix: 2 mL (200,000/mL x 11mL) 1.1x106 cells/mL	Add: 9 mL (11 mL - 2mL)

- 14. In a new sterile 50-mL conical tube, mix the calculated volumes of cell suspension and Basal Medium needed to obtain a final volume of 11 mL Astrocyte Seeding Suspension.
- 15. Add appropriate volumes of BrainFast Astrocyte (1:1000) and fetal bovine serum (1:100) to the Astrocyte Seeding Suspension per the Seeding Medium recipe.
 - a. Example: add 11 μ L of BrainFast Astro and 110 μ L fetal bovine serum to 11 mL of Astrocyte Seeding Suspension.
- 16. Mix completely and transfer 100 μ L/well (Ex: 20,000 cells/well) of the final Seeding Medium into a PDL-coated 96-well plate using a multi-channel pipette or liquid handler. Do not move or agitate the plate throughout the duration of the seeding process as this may lead to uneven attachment.
- 17. After seeding, let the plate rest for 10 minutes before moving it to allow the cells to settle to the bottom of the well. After 10 minutes, very gently transfer the plate to a humidified incubator at 37°C with 5% CO₂. The day of cell plating (today) is designated as Day 0.
 - *Note: The entire thawing and plating process should not exceed 1 hour. Post-thaw viability and overall cell health could be severely impacted and lead to an unsuccessful culture if the whole process is too long.

Day 1 Medium Replacement

- 18. On Day 1 (24 hours after seeding), prepare fresh Day 1 Medium.
- 19. Carefully remove all 100 μ L medium/well and gently replace with 100 μ L/well Day 1 Medium. Complete one row or column at a time to ensure the wells do not dry out during the medium replacement process.

Day 4 Medium Addition

- 20. On Day 4 (96 hours after seeding), prepare fresh Day 4 Medium.
- 21. Gently add 100 μ L/well of Day 4 medium to the entire plate for a final total volume of 200 μ L/well.

Day 7 and Onwards Medium Changes

- 22. Change half the medium (100 μ L/well) twice weekly using Basal Medium (made in Step 2).
 - a. Note: Addition of low concentration SK (0.1X-0.5X) in medium may be helpful for long-duration cultures; please contact support@brainxell.com for further assistance.
 - b. Gently remove 100 μ L/well and slowly add 100 μ L/well Basal Medium to the entire plate.

The astrocytes can be maintained viable and adherent in culture under the above conditions for at least 3 weeks post-seeding.

Astrocyte Media Compositions

Step 2		Component	Stock Conc.	Final	1 Plate	2 Plate	5 Plate
				Conc.	Volume	Volume	Volume
E	1	DMEM/F12 Medium	1 X	0.5X	24 mL	48 mL	120 mL
Medium	2	Neurobasal Medium	1 X	0.5X	24 mL	48 mL	120 mL
	3	N-2 Supplement	100X	1 X	500 μL	1 mL	2.5 mL
Basal	4	GlutaMAX	200 mM	0.5 mM	125 μL	250 μL	625 μL
Ğ	Basal Medium stock can be stored for up to 3 weeks at 4°C.						

Step 15		Component	Stock Conc.	Final Conc.	1 Plate Volume	2 Plate Volume	5 Plate Volume
Seeding Aedium 5	Astrocyte Seeding Suspension (Step 14)	1X	1X	11 mL	22 mL	55 mL	
	BrainFast Astro	1000X	1 X	11 μL	22 μL	55 μL	
0, 1	3	Fetal Bovine Serum	1X	1%	110 μL	220 μL	550 μL

Step 18		Component	Stock Conc.	Final Conc.	1 Plate Volume	2 Plate Volume	5 Plate Volume
- <u>E</u>	1	Basal Medium	1X	1X	11 mL	22 mL	55 mL
Day Mediu	2	BrainFast Astro	1000X	1 X	11 μL	22 μL	55 μL
_ >	3	BrainFast SK	1000X	1X	11 μL	22 μL	55 μL

Step 20		Component	Stock Conc.	Final Conc.	1 Plate Volume	2 Plate Volume	5 Plate Volume
ay 4 edium	1	Basal Medium	1X	1X	11 mL	22 mL	55 mL
Day	3	BrainFast Astro	1000X	1 X	11 μL	22 μL	55 μL