



---

### GENERAL GUIDELINES

#### RECOMMENDED FOR USE

The ComforTrac Cervical Traction Device is intended to treat a musculoskeletal or neurological impairment of the cervical spine to relieve pain, relax muscle spasms, and decompress spinal structures.

#### NOT RECOMMENDED FOR USE

Cervical Traction should be avoided in any condition of the cervical spine where movement can aggravate the condition or result in spinal instability, spinal injury and/or nerve root injury at risk for causing paralysis or ischemia.

#### SETUP:

1. Determine the patient's body weight. Tension should be up to 7% - 10% of body weight.
2. Determine the amount of neck flexion that is necessary and adjust the traction machine accordingly:  
Upper cervical region = 10 degrees (slide stand in upper most position)  
Middle cervical region = 15 degrees (slide stand in middle position)  
Lower cervical region = 20 degrees (slide stand in lower most position)
3. Remove any jewelry (especially earrings), eyeglasses, and anything else in the cervical region that may get in the way or create discomfort for the patient. Loosen the shirt collar or have the patient change into a gown in order to better expose the neck region. (Cervical traction is only effective if applied if the traction unit is in direct contact with the skin of the neck, not over clothing.)
4. Instruct the patient to lay supine on the treatment table with the neck between the neck wedges of the cervical traction device. Adjust the neck wedges appropriately so that they are snug around the patient's neck. If necessary have the patient sit up and adjust the neck wedges and then have the patient lie back down.
5. Use the hand pump to increase pressure to appropriate level. (Please refer to the protocols within this section for guidelines on the tension amount and the treatment time for specific diagnoses.)



---

### ACUTE CERVICAL PAIN

Cervical traction may be used for acute cervical pain as long as the cervical spine has been cleared for any of the Recommended for Use information mentioned within the General Guidelines section. Traction may also be more beneficial when used with other modalities such as cold, heat, ultrasound, etc. to reduce pain. Initially, apply tension only up to patient tolerance as aggressive stretching of the cervical spine may create muscular spasm and pain in this condition.

#### SUGGESTED SETTINGS

##### GENERAL:

Tension: 10 - 20 pounds\*

Tension type: Sustained tension

Treatment Time: 10-15 minutes

##### VARIATIONS:

Upper cervical region = 10 degrees flexion (slide stand in upper most position)

Middle cervical region = 15 degrees flexion (slide stand in middle position)

Lower cervical region = 20 degrees flexion (slide stand in lower most position)

Intermittent tension setting: 2:1 up to 4:1 ratio of tension to rest (e.g. 2-4 minutes on, 1 minute rest)

**SESSIONS PER DAY:** 1-3

---

### CAPSULITIS OF CERVICAL SPINE JOINTS

Cervical traction may be effective in helping to treat Capsulitis of the cervical spine. Be sure that the cervical spine has been cleared for any of the Recommended for Use information within the General Guidelines section. Traction may also be more beneficial when used with other modalities such as cold, heat, ultrasound, etc.

#### SUGGESTED SETTINGS

##### GENERAL:

Tension: 10 - 30 pounds (may be increased up to 7-10% of patient body weight)

Tension type: Intermittent tension; 2:1 up to 4:1 ratio of tension to rest (e.g. 2-4 minutes on, 1 minute rest)

Treatment Time: 10-15 minutes

##### VARIATIONS:

Upper cervical region = 10 degrees flexion (slide stand in upper most position)

Middle cervical region = 15 degrees flexion (slide stand in middle position)

Lower cervical region = 20 degrees flexion (slide stand in lower most position)

Sustained tension setting: Apply sustained tension to cervical spine with periodic rests during the 10-15 minute treatment time

**SESSIONS PER DAY:** 1-3



---

### CERVICAL SPONDYLOSIS

Cervical traction may be used for Radiculopathy as long as the cervical spine has been cleared for any of the Recommended for Use information within the General Guidelines section. Traction may also be more beneficial when used with other modalities such as cold, heat, ultrasound, etc.

#### SUGGESTED SETTINGS

##### GENERAL:

Tension: 10 - 30 pounds (may be increased up to 7-10% of patient body weight)

Tension type: Sustained tension

Treatment Time: 10-15 minutes

##### VARIATIONS:

Upper cervical region = 10 degrees flexion (slide stand in upper most position)

Middle cervical region = 15 degrees flexion (slide stand in middle position)

Lower cervical region = 20 degrees flexion (slide stand in lower most position)

Intermittent tension setting: 2:1 up to 4:1 ratio of tension to rest (e.g. 2-4 minutes on, 1 minute rest)

**SESSIONS PER DAY:** 1-3

---

### CHRONIC CERVICAL PAIN

Cervical traction may be used for chronic cervical pain as long as the cervical spine has been cleared for the Recommended for Use information within the General Guidelines section. Traction may also be more beneficial when used with other modalities such as cold, heat, ultrasound, etc.

#### SUGGESTED SETTINGS

##### GENERAL:

Tension: 10 - 30 pounds (may be increased up to 7-10% of patient body weight)

Tension type: Intermittent tension; 2:1 up to 4:1 ratio of tension to rest (e.g. 2-4 minutes on, 1 minute rest)

Treatment Time: 10-15 minutes

##### VARIATIONS:

Upper cervical region = 10 degrees flexion (slide stand in upper most position)

Middle cervical region = 15 degrees flexion (slide stand in middle position)

Lower cervical region = 20 degrees flexion (slide stand in lower most position)

Sustained tension setting: Apply sustained tension to cervical spine with periodic rests during the 10-15 minute treatment time

**SESSIONS PER DAY:** 1-3



---

### DEGENERATIVE DISC DISEASE (DDD)

Degenerative disc disease (DDD) often leads to radiculopathy and cervical pain. Cervical traction alleviates pressure on the diseased discs to promote decreased pressure on peripheral nerve roots and enhance nutrient delivery to the disc. Be sure that the cervical spine has been cleared for any of the Recommended for Use information within the General Guidelines section. Traction may also be more beneficial when used with other modalities such as cold, heat, ultrasound, etc.

#### SUGGESTED SETTINGS

##### GENERAL:

Tension: 10 - 30 pounds (may be increased up to 7-10% of patient body weight)

Tension type: Sustained tension Treatment Time: 10-15 minutes

##### VARIATIONS:

Upper cervical region = 10 degrees flexion (slide stand in upper most position)

Middle cervical region = 15 degrees flexion (slide stand in middle position)

Lower cervical region = 20 degrees flexion (slide stand in lower most position)

Intermittent tension setting: 2:1 up to 4:1 ratio of tension to rest (e.g. 2-4 minutes on, 1 minute rest)

**SESSIONS PER DAY:** 1-3

---

### FACET JOINT DYSFUNCTION

Cervical traction may be effective to increase the mobility of the cervical facet joints as long as the cervical spine has been cleared for any of the Recommended for Use information within the General Guidelines section. Traction may also be more beneficial when used with other modalities such as cold, heat, ultrasound, etc.

#### SUGGESTED SETTINGS

##### GENERAL:

Tension: 10 - 30 pounds (may be increased up to 7-10% of patient body weight)

Tension type: Intermittent tension; 2:1 up to 4:1 ratio of tension to rest  
(e.g. 1-2 minutes on, 30 second rest)

Treatment Time: 10-15 minutes

##### VARIATIONS:

Upper cervical region = 10 degrees flexion (slide stand in upper most position)

Middle cervical region = 15 degrees flexion (slide stand in middle position)

Lower cervical region = 20 degrees flexion (slide stand in lower most position)

Sustained tension setting: Apply sustained tension to cervical spine with periodic rests during the 10-15 minute treatment time

**SESSIONS PER DAY:** 1-3



---

### HEADACHES

Cervical traction may be used for Radiculopathy as long as the cervical spine has been cleared for any Recommended for Use information within the General Guidelines section. Traction may also be more beneficial when used with other modalities such as cold, heat, ultrasound, etc.

#### **SUGGESTED SETTINGS**

##### **GENERAL:**

Tension: 10 - 30 pounds (may be increased up to 7-10% of patient body weight)

Tension type: Intermittent tension; 2:1 up to 4:1 ratio of tension to rest (e.g. 2-4 minutes on, 1 minute rest)

Treatment Time: 10-15 minutes

##### **VARIATIONS:**

Upper cervical region = 10 degrees flexion (slide stand in upper most position)

Middle cervical region = 15 degrees flexion (slide stand in middle position)

Lower cervical region = 20 degrees flexion (slide stand in lower most position)

Sustained tension setting: Apply sustained tension to cervical spine with periodic rests during the 10-15 minute treatment time

**SESSIONS PER DAY:** 1-3

---

### HERNIATED CERVICAL DISC

A herniated cervical disc often leads to radiculopathy and cervical pain. Cervical traction alleviates pressure on the diseased discs to promote decreased pressure on peripheral nerve roots and enhance nutrient delivery to the disc. Be sure that the cervical spine has been cleared for any of the Recommended for Use information within the General Guidelines section. Traction may also be more beneficial when used with other modalities such as cold, heat, ultrasound, etc.

#### **SUGGESTED SETTINGS**

##### **GENERAL:**

Tension: 10 - 30 pounds (may be increased up to 7-10% of patient body weight)

Tension type: Sustained tension

Treatment Time: 10-15 minutes

##### **VARIATIONS:**

Upper cervical region = 10 degrees flexion (slide stand in upper most position)

Middle cervical region = 15 degrees flexion (slide stand in middle position)

Lower cervical region = 20 degrees flexion (slide stand in lower most position)

Intermittent tension setting: 2:1 up to 4:1 ratio of tension to rest (e.g. 2-4 minutes on, 1 minute rest)

**SESSIONS PER DAY:** 1-3



---

### JOINT HYPOMOBILITY — CERVICAL SPINE

Cervical traction increases joint mobility of the cervical spine for patients with hypomobility of the cervical spine. Be sure that the cervical spine has been cleared for any of the Recommended for Use information within the General Guidelines section. Traction may also be more beneficial when used with other modalities such as cold, heat, ultrasound, etc.

#### SUGGESTED SETTINGS

##### GENERAL:

Tension: 10 - 30 pounds (may be increased up to 7-10% of patient body weight)

Tension type: Intermittent tension; 2:1 up to 4:1 ratio of tension to rest (e.g. 2-4 minutes on, 1 minute rest)

Treatment Time: 10-15 minutes

##### VARIATIONS:

Upper cervical region = 10 degrees flexion (slide stand in upper most position)

Middle cervical region = 15 degrees flexion (slide stand in middle position)

Lower cervical region = 20 degrees flexion (slide stand in lower most position)

Sustained tension setting: Apply sustained tension to cervical spine with periodic rests during the 10-15 minute treatment time

**SESSIONS PER DAY:** 1-3

---

### MUSCULAR SPASM OF CERVICAL MUSCLES

Cervical traction may be effective in helping to relax muscular spasms of the cervical spine. Be sure that the cervical spine has been cleared for any of the Recommended for Use information within the General Guidelines section. Traction may also be more beneficial when used with other modalities such as cold, heat, ultrasound, etc.

#### SUGGESTED SETTINGS

##### GENERAL:

Tension: 10 - 30 pounds (may be increased up to 7-10% of patient body weight)

Tension type: Sustained tension

Treatment Time: 10-15 minutes

##### VARIATIONS:

Upper cervical region = 10 degrees flexion (slide stand in upper most position)

Middle cervical region = 15 degrees flexion (slide stand in middle position)

Lower cervical region = 20 degrees flexion (slide stand in lower most position)

Intermittent tension setting: 2:1 up to 4:1 ratio of tension to rest (e.g. 2-4 minutes on, 1 minute rest)

**SESSIONS PER DAY:** 1-3



---

### OSTEOARTHRITIS — CERVICAL SPINE

Cervical traction may benefit patients with osteoarthritis by promoting increased movement of the cervical spine and decreased pressure on the facet joints. Be sure that the cervical spine has been cleared for any of the Recommended for Use information within the General Guidelines section. Traction may also be more beneficial when used with other modalities such as cold, heat, ultrasound, etc.

#### **SUGGESTED SETTINGS**

##### **GENERAL:**

Tension: 10 - 30 pounds (may be increased up to 7-10% of patient body weight)

Treatment Time: 10-15 minutes

##### **VARIATIONS:**

Upper cervical region = 10 degrees flexion (slide stand in upper most position)

Middle cervical region = 15 degrees flexion (slide stand in middle position)

Lower cervical region = 20 degrees flexion (slide stand in lower most position)

Intermittent tension setting: 2:1 up to 4:1 ratio of tension to rest (e.g. 2-4 minutes on, 1 minute rest)

**SESSIONS PER DAY:** 1-3

---

### RADICULOPATHY

Cervical traction effectively reduces radiculopathy resulting from impingement of cervical nerve roots. Be sure that the cervical spine has been cleared for any of the Recommended for Use information within the General Guidelines section. Traction may also be more beneficial when used with other modalities such as cold, heat, ultrasound, etc.

#### **SUGGESTED SETTINGS**

##### **GENERAL:**

Tension: 10 - 30 pounds (may be increased up to 7-10% of patient body weight)

Tension type: Sustained tension

Treatment Time: 10-15 minutes

##### **VARIATIONS:**

Upper cervical region = 10 degrees flexion (slide stand in upper most position)

Middle cervical region = 15 degrees flexion (slide stand in middle position)

Lower cervical region = 20 degrees flexion (slide stand in lower most position)

Intermittent tension setting: 2:1 up to 4:1 ratio of tension to rest (e.g. 2-4 minutes on, 1 minute rest)

**SESSIONS PER DAY:** 1-3





## Cervical Traction Protocols



---

### WHIPLASH

Cervical traction may be an effective for Whiplash as long as the cervical spine has been cleared for any of the Recommended for Use information within the General Guidelines section. Traction may also be more beneficial when used with other modalities such as cold, heat, ultrasound, etc.

#### SUGGESTED SETTINGS

##### GENERAL:

Tension: 10 - 30 pounds (may be increased up to 7-10% of patient body weight)

Tension type: Sustained tension

Treatment Time: 10-15 minutes

##### VARIATIONS:

Upper cervical region = 10 degrees flexion (slide stand in upper most position)

Middle cervical region = 15 degrees flexion (slide stand in middle position)

Lower cervical region = 20 degrees flexion (slide stand in lower most position)

Intermittent tension setting: 2:1 up to 4:1 ratio of tension to rest (e.g. 2-4 minutes on, 1 minute rest)

#### SESSIONS PER DAY: 1-3

---



Distributed by:

**RS Medical**<sup>®</sup>  
*Pain & Function Solutions<sup>SM</sup>*

ComforTrac Inc.  
3930 Horseshoe Bend Road  
Keysville, VA 23947

866.378.4672  
703.891.0455  
703.288.3304 fax

info@comfortrac.net  
[www.comfortrac.net](http://www.comfortrac.net)

14001 SE First Street  
Vancouver, WA 98684  
[www.rsmedical.com](http://www.rsmedical.com)

360.892.0339  
800.683.0353  
360.896.2566 fax