

Tunstall

Liberty³⁰⁰ Two-Way Vibrating Pendant User Guide



Product overview

The Liberty³⁰⁰ Two-Way Vibrating Pendant is our new frequency hopping spread spectrum pendant with updated receiver integration.

Your pendant (personal radio trigger) is wirelessly connected to your Liberty³⁰⁰ home unit and enables you to generate a call for help from anywhere in your home or garden.

The pendant is bi-directional, meaning that when it sends a signal to the Liberty³⁰⁰ home unit and the signal is acknowledged by the alarm, the pendant will vibrate and the red LED will also illuminate.

The pendant should be worn at all times and can be worn around the neck or on the wrist using the attachments supplied.

Your pendant is water resistant (to IP67 standard) up to one metre for 30 minutes. It can be worn in the shower however wearing it in the bath should be avoided where possible.

Programming to the Liberty³⁰⁰

Once received, the Liberty³⁰⁰ Two-Way Vibrating Pendant will require programming to your alarm, by following the steps below.

1. Call **1800 046 981** and advise the Operator you have a new pendant to be programmed.
2. Give the Operator your contact number and name for confirmation.
3. Confirm the Serial Number of the pendant located underneath the barcode on the back of your pendant.



The Operator will then send information through to your alarm, once completed, you will be asked to test the pendant by pressing the button.

Your new Liberty³⁰⁰ Two-Way Vibrating Pendant is now active and ready to use.

Cleaning

If required, the Liberty³⁰⁰ Two-Way Vibrating Pendant may be cleaned periodically with a damp cloth and mild detergent or an alcohol wipe.

Technical details

Dimensions (W x H x D)	38mm x 40mm x 12mm
Weight	25g
Battery type	Lithium
Battery life	10 year lithium quiescent standby
Radio frequency	915-928MHz frequency range
Radio range	100-300m in open air (subject to environmental conditions)
Waterproof rating	IP67 / Up to 1m for 30mins
Operating temp range	Ambient temperature range 0°C to +50°C
Radio (main unit and pendant):	AS/NZS 4268.2012
Design, Manufacture:	ISO9001:2000

Manufactured in Australia

tunstallhealthcare.com.au

Our policy of continual development means that product specifications and appearance may change without notice. Tunstall does not accept any responsibility for any errors and omissions contained within this document.

© 2017 Tunstall Group Ltd.

®Tunstall is a registered trade marks

Tunstall Australasia Pty Ltd
1/56 Lavarack Ave, Eagle Farm QLD 4009
ABN 44 059 121 863
ACN 059 121 863

Tunstall New Zealand Ltd
2/65 Chapel St, Tauranga 3110
Business No. 3502431

Tunstall

Life Changing, Life Saving

#8079 v1.0