



PLYMOUTH
MEDICAL

DERMAL REJUVENATION PROCEDURE CHECKLIST

Equipment purchasable from Plymouth Medical

- o (1) EmCyt^e Centrifuge
- o (1) Pure PRP II kit (GS30 PURE II 4mL PRP injectate, GS 60 Pure II PRP 7mL PRP injectate)
- o BLT in lipoderm 20%/6%/4% 100g Jar MU
- o Tulip Cannulas
- o Nitronox machine
- o DP DermaceuticalsTM Vitamin Rich RepairTM (30ml)
- o DP DermaceuticalsTM Antioxidant CocktailTM (30ml)
- o DP DermaceuticalsTM Tri-Phase CleanserTM (30ml)
- o DP DermaceuticalsTM 3D Sculptured Mask

Kit	Description
GS30-PURE-II PurePRP [®] II Platelet Concentrating System 30mL	• 3- 30mL Syringe • 1- 12mL Syringe • 1- 3mL Syringe with 18G x 1" Needle • 1- Concentrating Device 30mL • 1- Concentrating Accessory 30mL • 1- Blood Draw Accessory • 1- Sterile Drape • 1- Citrate Anticoagulant 20mL • 1-18G x 1" Needle • 4- Alcohol Pre-Pads • 6- Panel Labels
GS60-PURE-II PurePRP [®] II Platelet Concentrating System 60mL	• 2- 60mL Syringe • 1- 30mL Syringe • 1- 12mL Syringe • 1- 3mL Syringe with 18G x 1" Needle • 1- Concentrating Device 60mL • 1- Concentrating Accessory 60mL • 1- Blood Draw Accessory • 1- Sterile Drape • 1- Citrate Anticoagulant 20mL • 1-18G x 1" Needle • 4- Alcohol Pre-Pads • 6- Panel Labels

PLYMOUTH MEDICAL

135 PLYMOUTH ST, SUITE 410, BROOKLYN, NY, 11201 • TEL 888-392-5076 • FAX 888-392-1455 Federal ID# 45-3973550
WWW.PLYMOUTHMEDICAL.COM



Required supplementary items

1. Preparation & procedure supplies (additional items to purchase)

- o (1) Pair Sterile Gloves (per physician/assistant) SU
- o (1) Sterile Surgical Skin Marker SU
- o (1) Bacteriostatic Saline vial MU
- o (5) Sterile Gauze Pad 4x4 SU
- o 1mL Syringes SU (GENERAL RULE OF THUMB: 1mL luer syringes for each ML of PRP created for injections)

2. Needle recommendations/options

- o Tulip Gauge Injection Cannula (for Tear trough and Lip injection)
- o 27 or 30-gauge needle for Dermal injections

3. Anesthetics/Analgesics recommendation

- o BLT in lipoderm 20%/6%/4%
- o Nitronox Device
- o Lidocaine 1% (Lido Block)

Processing Video for No Burn Protocol (Micro needling with Injections / Hair / P shot/ O Shot)

For GS30 Pure II PRP kits and GS60 Pure II PRP kits:

<https://www.youtube.com/watch?pbjreload=10&v=lo9waUiTkBA>

PLYMOUTH MEDICAL

135 PLYMOUTH ST, SUITE 410, BROOKLYN, NY, 11201 • TEL 888-392-5076 • FAX 888-392-1455 Federal ID# 45-3973550
WWW.PLYMOUTHMEDICAL.COM