

# Nitronox Plus® 0-70%, Nitronox Plus® 50/50 & Pro-Nox™ Comparison



PLYMOUTH  
MEDICAL

Features	Nitronox Plus® 0-70%	Nitronox Plus® 50/50	Pro-Nox™
N <sub>2</sub> O can be adjusted from 0-70%	✓		
Can deliver 100% O <sub>2</sub> - flushes patient's system of any residual nitrous oxide to help patients recover faster	✓	✓	
Visual Gas Gauges - real time confirmation of gas supply and mixture pressure	✓	✓	
Key Lock System - provides security and prevents device from being used without the key	✓	✓	
Function Test Button - easy access to perform periodic safety and function tests	✓	✓	
2-Cylinder & 4-Cylinder Options Available	✓	✓	
Internal, But Easily Accessible Demand Valve For Self Maintenance	✓	✓	
3 Year Warranty (Headpiece)	✓	✓	
Manufactured & Serviced in the USA	✓	✓	
Patient Self-Administered	✓	✓	✓
Audible Low Gas Pressure Alarm - alerts staff when either O <sub>2</sub> or N <sub>2</sub> O gas supply pressure is low	✓	✓	✓
Emergency Air Intake - allows patient to pull in room air through breathing circuit should gas supply run out	✓	✓	✓
Delivers 50% Mixture of N <sub>2</sub> O & O <sub>2</sub>	✓	✓	✓

Nitronox Plus® 0-70%



Nitronox Plus® 50/50




Pro-Nox™



Headpiece Exterior




Headpiece Interior



Headpiece Exterior



Headpiece Interior



Headpiece Exterior



Headpiece Interior

