

# Microbiology

## Germ-Free Rodent Testing (Call for Details)

Service:	Description:	Sample Type:
GF-1	Inoculation of three permissive culture media (2 solid media, 1 broth media), incubated aerobically. For the detection of aerobic bacterial growth. Includes identification of growth by MALDI-TOF MS*	Feces or Swabs
GF-2	Inoculation of three permissive culture media (2 solid media, 1 broth media), incubated aerobically. Plus, inoculation of 1 solid media anaerobically for 3 days. For the detection of aerobic and anaerobic bacterial growth. Includes identification of growth by MALDI-TOF MS*	Feces or Swabs
Axenic Animal Sterility w/ Anaerobic Culture	Inoculation of 8 permissive culture media (4 solid media, 4 broth media), incubated aerobically. Plus 2 solid media incubated anaerobically for 10 days. For the detection of aerobic bacteria, anaerobic bacteria, fastidious bacteria, filamentous fungi and yeasts. Includes identification of growth by MALDI-TOF MS*	Feces or Swabs