



KVP Rehabilitation

Balto® Flex

A Brace for Neurological Disorders



The Balto® Flex is used to help dogs with neurological disorders that cause severe proprioceptive deficits, in which dogs are unable to sense the position of their limb on the ground. The brace can also be used in cases of sciatic nerve damage to ensure the neurological problem does not lead to orthopedic damage.

ITEM #	SIZE	WEIGHT
BTPULLXS	XS	10 - 20 lbs
BTPULLS	S	20 - 45 lbs
BTPULLM	M	45 - 85 lbs
BTPULLL	L	over 85 lbs