

Disability & Hearing Impairment Solutions

Bellman Visit helps people with disability and impaired hearing to detect important signals and alarms in their home, like for instance the doorbell, telephone, crying baby, smoke alarm etc.

Bellman Visit was developed to bring attention and warn people of events in their natural environments. It can be interfaced with Tunstall alarms and sensors to provide additional support to connected care users, or a unique alerting solution for onsite carers.



Bellman Visit Door Transmitter

The Visit Door Transmitter helps you notice when the doorbell rings, even if you are several rooms away. It signals the Visit receiver that notifies you with your own choice of flashes, sounds and vibrations.



Bellman Visit Baby Cry Transmitter

The Baby Cry transmitter is activated by crying or screaming when placed near the baby. Used in conjunction with Tunstall alarms, it may also be used for other home care applications.



Bellman Visit Smoke Alarm Transmitter

The Visit smoke alarm transmitter helps the user detect smoke and smoldering fire at an early stage. It is mounted on the ceiling and sends a radio signal to the Visit receiver for alerting once activated by smoke.



Bellman Visit Telephone Transmitter

The Visit telephone transmitter helps you to notice when the telephone rings. When activated the Visit receiver notifies you with your own choice of flashes, sounds and vibrations.

It can also be customised for many other applications and Tunstall sensors.



Bellman Visit Push Button Transmitter

The Visit push button transmitter works as a wireless doorbell and can be used as a personal alarm notifying the other residents of the comings and goings in your home or if the wearer needs assistance.



Bellman Door Contact Mat

The contact mat together with the Visit multi-purpose transmitter is used to notify when someone gets out of bed or enters/leaves a room. When a person steps on it, the Visit transmitter signals the Visit receivers that notifies you. The doormat can also interface with a Tunstall alarm and be used as an activity sensor.



Bellman Visit Portable Receiver

The Visit portable receiver alerts the user with a range of different lights and sounds when various Visit transmitters send signals. When interfaced with Tunstall sensors the portable receiver becomes a convenient alerting device for onsite carers.



Bellman Alarm Clock

The alarm clock awakes you with a loud sound or flash lights and bed shaker vibrations when its time to wake up or when a Visit transmitter is activated. In addition Tunstall sensors can also interface with the alarm clock to wake up sleeping carers.



Bellman Visit Pager Receiver

The Visit pager notifies the user with vibrations and lighted signals when the Visit transmitter is activated. It is small, light-weight, can be carried in the pocket or on the belt. The Tunstall sensor can be interfaced to trigger the pager alerting the onsite carer. During the night, you can connect a bed shaker to the charger.



Bellman Bed Shaker

The bed shaker is placed under the pillow or mattress to alert a heavy sleeper of calls and alarms. It works with all Visit receivers and requires no internal battery.



Bellman Visit Flash Receiver

The Visit flash receiver attracts the users attention with bright flashes with different colored LED lights when Visit transmitter sends signals. This easy to use portable flash receiver comes with a mounted table stand and a wall bracket accessory.



Tunstall Australasia Pty Ltd

ABN 44 059 121 863

Locked Bag 1, 985 Kingsford Smith Dr
Eagle Farm, QLD 4009

Tel: 1800 708 758

Email: customer.service@tunstallhealthcare.com.au

www.tunstallhealthcare.com.au

Tunstall New Zealand Ltd

Business No. 3502431

PO Box 13153
Tauranga 3110

Tel: 0800 488 678

Email: info@tunstall.co.nz

www.tunstall.co.nz



Bellman & Symfon®

Tunstall