

TF Thomas Henry Fonzi
Whidbey Island Clinic - Oak Harbor

Detectors
Operators
Providers
Facilities
Settings

Test Results
Patients
Start New Test
Scan Test Result

⚙️

Settings

Organization Details

Patient Portal & Report Letters

Use these settings to configure your organization's Patient Portal and Report Letters. The Patient Portal must be On to enable Report Letters.

Enable Patient Portal On Off

The Patient Portal allows your organization's patients to view and print their test results from a secure web portal.

Invite New Patients On Off

Automatically invite new patients to create a Patient Portal account.

Report Letters

Report Letters provide your patients with a formal copy of their COVID-19 test results.

Report Letter Preview

Options

Report Letter Template

ReportLetterSansSerif.html (download sample PDF or download template)

ReportLetterSerif.html (download sample PDF or download template)

Upload New Template

Official Signature

Letters must include the signature of an official representative of your organization. Please enter the name and title of the signee for your organization here to generate the signature.

First Name

Last Name

Credentials

Title

Thomas Henry Fonzi

Thomas Henry Fonzi, NPP
NP

Confirm Signature

Introducing our Exclusive Detector Portal



The AscencioDx[®] Molecular Detector is a rapid, highly accurate, easy-to-use and affordable clinical diagnostic platform for performing molecular analysis of patient samples in point-of-care (POC) settings using The AscencioDx COVID-19 Test.

The AscencioDx[®] COVID-19 Test and The AscencioDx[®] Molecular Detector has not been FDA cleared or approved, but has been authorized for emergency use by FDA under an EUA for use by authorized laboratories; This product has been authorized only for the detection of nucleic acid from SARS-CoV-2, not for any other viruses or pathogens; and The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug, and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.



[Learn more >>](#)

[Anavasi Diagnostics >>](#)

Begin using your AscencioDx[®] Molecular Detector

An Anavasi Diagnostics account allows you to manage the patients and operators who use your detector, review the testing history of patients, and simplifies the day-to-day use of your detector.

[Login](#)

[Create an Account](#)

Anavasi Diagnostics fully encrypts your data and securely sends it to relevant federal, state, and local health authorities via [CDC ReportStream](#).

- Comprehensive report generation and customizable "Passport" style report templates with option to upload your own template

- Results are securely transmitted from the detector to our cloud-based portal via a QR-Code
- No Bluetooth required or app to download
- HIPAA compliant
- Automatically disseminates results to the CDC ReportStream, which in turn reports to in-program State and Local Health Agencies
- Eliminate the hassle of paper or cumbersome spreadsheets while easily and securely managing patient records on your dashboard

Note: This product has not been FDA cleared or approved, but has been authorized for emergency use by FDA under an EUA for use by authorized laboratories. This product has been authorized only for the detection of nucleic acid from SARS-CoV-2, not for any other viruses or pathogens. The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug, and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

Anavasi Diagnostics
Redmond, WA USA
www.anavasidx.com
info@anavasidx.com
PH: 800-262-8274



TF Thomas Henry Fonzi
Whidbey Island Clinic - Oak Harbor

Detectors Operators Providers Facilities Settings

Test Results
Patients
Start New Test
Scan Test Result

Settings

Organization Details

Patient Portal & Report Letters

Your account contact@whidbeyclinic.com gives you administrative privileges and permissions for your organization's information on the Anavasi Diagnostics website. As an administrator, you have the ability to run tests and to manage patients, detectors, test operators, providers, and facilities on behalf of your organization. These settings apply for all administrators, operators, and patients in your organization. For settings specific to your account click on your username above.

Organization Name

Whidbey Island Clinic Save Name

Show Quality Controls Show Hide
Controls whether Quality Controls (QCs) may be assigned to test results carried out at Sunshine State facilities.

Special Features ACTIVE

These features let you generate synthetic test results and detectors for training or demonstration purposes. Synthetic test results are never reported to public health authorities.

Mark all new test result QR code scans as synthetic ACTIVE

Enable generating synthetic test results from Detectors page ACTIVE

Enable adding synthetic detectors from Add New Detector ACTIVE

Request Assistance INACTIVE

Enabling this feature will send a request to Anavasi to aid you in configuring your portal. Doing this will grant an Anavasi representative the ability to log in and make modifications to your portal.

Logo

Upload your organization's logo. This logo will be used on your organization's Patient Portal and in result Report Letters

Select File

Organization Contact Details

Provide your organization's main address, telephone number, and fax number.

Address
32550 State Route 20

City
Miami

State
Washington

Zip Code
98277

Phone Number
(360) 555-5555

FAX Number
(360) 555-7777


Email
contact@whidbeyclinic.com

Timezone
US/Pacific

Web Address
http://www.WhidbeyClinic.com

Save Contact Details

- Value added portal exclusively available for The AscencioDx Molecular Detector at no additional cost
- Set up your organization and manage detectors, locations, operators, physicians and patient data on your dashboard



TF Thomas Henry Fonzi
Whidbey Island Clinic - Oak Harbor

Detectors Operators Providers Facilities Settings

Test Results
Patients
Start New Test
Scan Test Result

Test Results – Whidbey Island Clinic

Below are the test results from all detectors in the Whidbey Island Clinic - Oak Harbor facility. New results are added when test result QR codes are scanned.

Export Show Assigned patients Unassigned patients 07/06/2023 - 10/04/2023

PATIENT	TEST RESULT	DATE / DETECTOR	PROVIDER	PATIENT PORTAL
<p>Ratner, Kim DOB: 01 / 02 / 1988 @ . - () -</p>	<p>☹️ POSITIVE COVID-19 Virus Detected</p>	<p>06/20/23 12:28 PM SN 5668986/1 Oak Harbor</p>	<p>Tim Barton, NP</p>	<p>Regenerate Letter Email / Invite Print / Download</p>
<p>Ratner, Kim DOB: 01 / 02 / 1988 @ . - () -</p>	<p>😊 NEGATIVE COVID-19 Virus Not Detected</p>	<p>06/20/23 12:34 PM SN 5668986/2 Oak Harbor</p>	<p>Denise Marie Andrews, Supervisor</p>	<p>Generate Letter Email / Invite Print / Download</p>
<p>Ratner, Kim DOB: 01 / 02 / 1988 @ . - () -</p>	<p>☹️ POSITIVE COVID-19 Virus Detected</p>	<p>06/20/23 12:42 PM SN 5668986/3 Oak Harbor</p>	<p>Keith Kelly, NP</p>	<p>Generate Letter Email / Invite Print / Download</p>
<p>Adcock, Marilyn DOB: 04 / 04 / 2004 @ . - () -</p>	<p>😊 NEGATIVE COVID-19 Virus Not Detected</p>	<p>06/22/23 2:00 PM SN 5668986/4 Oak Harbor</p>	<p>Tanya Lloyd, NP</p>	<p>Generate Letter Email / Invite Print / Download</p>
<p>Farrow, Xaria DOB: 01 / 01 / 2001 @ . - () -</p>	<p>☹️ POSITIVE COVID-19 Virus Detected</p>	<p>06/22/23 2:06 PM SN 5668986/5 Oak Harbor</p>	<p>Conrad Richards, NP</p>	<p>Generate Letter Email / Invite Print / Download</p>
<p>Farrow, Xaria DOB: 01 / 01 / 2001 @ . - () -</p>	<p>😊 NEGATIVE COVID-19 Virus Not Detected</p>	<p>06/24/23 12:55 PM SN 5668986/6 Oak Harbor</p>	<p>Tanya Lloyd, NP</p>	<p>Generate Letter Email / Invite Print / Download</p>
<p>Ratner, Kim DOB: 01 / 02 / 1988 @ . - () -</p>	<p>☹️ POSITIVE COVID-19 Virus Detected</p>	<p>06/24/23 12:58 PM SN 5668986/7 Oak Harbor</p>	<p>Tanya Lloyd, NP</p>	<p>Generate Letter Email / Invite Print / Download</p>

- A click of a button securely emails test results to a patient
- Provide individual portal access to patients so they have a secure record of all their COVID test results
- Results are also exportable in an excel format
- Patient information can be auto-populated by scanning their Driver's License -- reducing data entry errors and saving time.