

FrontLine™

FrontLine™ Powder-Free Nitrile Exam Gloves were developed to offer the following:

- Human skin pH 5.5 balanced coating for reduced risk of irritation and perspiration
- Advanced non-latex nitrile polymer offers superior stretch and comfort
- Fully textured surface for improved grip in all conditions
- UL certified NFPA 1999-2018 edition
- Exceeds Puncture Resistance of NFPA 1999-2018 edition
- **Fentanyl Permeation Testing:** FrontLine™ has zero breakthrough with Fentanyl Citrate 100mcg/2ml for 240 minutes (test limit) per ASTM D6978. This assures EMS, Law Enforcement, and Corrections maximum protection in the event of inadvertent contact with synthetic opiates. Test results on file.
- Our **EaseOff™** non-stick exterior surface prevents tapes including Duct Tape and even **Veni-Gard™** from sticking to the glove
- AQL of 0.65 exceeding industry standards
- **ExactFit™** glove formers for reduced hand fatigue & improved dexterity
- ISO 13485 certified by an accredited registrar.



Ordering information

Size	Catalog no.	Packaging
Small	FL9-1	100 gloves / box
Medium	FL9-2	100 gloves / box
Large	FL9-3	100 gloves / box
XLarge	FL9-4	100 gloves / box
XXLarge	FL9-5	90 gloves / box
XXXLarge	FL9-6	90 gloves / box

Size	Catalog no.	Packaging
Small	FL12-1	100 gloves / box
Medium	FL12-2	100 gloves / box
Large	FL12-3	100 gloves / box
XLarge	FL12-4	100 gloves / box
XXLarge	FL12-5	90 gloves / box
XXXLarge	FL12-6	90 gloves / box

All packaging – 10 boxes / case

Product specification

Specifications	ASTM D6319	NFPA 2018 Ed	FrontLine™	
			FL9	FL12
Length (in)	9.05"	n/a	9.5"	12"
Fingertip (mils)	1.9mil	n/a	6mil	6mil
Palm (mils)	1.9mil	n/a	5mil	5mil
Cuff (mils)	n/a	n/a	4mil	4mil
Tensile strength				
Before aging (MPa)	14	14	20	20
After aging (MPa)	14	14	18	18
Elasticity / elongation				
Before aging (%)	500%	500%	600	600
After aging (%)	400%	500%	510%	510%
Puncture resistance	n/a	1.01 lbs	1.9 lbs	1.9 lbs

