

HOSPITAL

Sanitize Your C-Com[®] Headsets & Other Devices at the Speed of Light



An average of 1 in 4 hospital devices are contaminated with pathogenic bacteria.

We can help fix that.



M00033A

A Safe & Easy Solution to Reduce the Spread of Bacteria



Fast and simple. No training required.




No chemicals and no damage to sensitive touchscreens.



Touch-free device removal that enables proper hand hygiene.



Proven to inactivate 99.9998% of MRSA without pre-cleaning.



www.carrotmedical.com | 425-318-8089



Break the Chain of Infection



Simple, fast, and effective device sanitization station



**Designed
for Hospitals.
Trusted
Everywhere.**

2600+
Deployments

Zero
Damage to
Touchscreens

350+
Hospitals Across
USA & CAN

720
Devices per
Hour Capacity

Efficacy Rate @ 20 Secs

- S. aureus (MRSA) - 99.99981%
- C. difficile spores - 99.97%
- S. enterica (Salmonella) - 99.999%
- E. faecalis (VRE) - 99.9992%
- E. coli - 99.9991%
- SARS-CoV-2 - 99.995%

Specifications

- **Size:** 26.7" x 16.5" x 12.3"
- **Tray:** 9.4" x 12.87" x 3.45"
- **Desktop Footprint:** 15.7" x 16.45"



Trusted by Over 350 Hospitals
to Sanitize Tablets, Headsets & More.



info@carrotmedical.com
www.carrotmedical.com
+1 425 318 8089